

.

---

**From:** michelle\_tomba@agilent.com  
**Sent:** Monday, July 26, 2010 6:10 PM  
**To:** Hanchett, James (DPH)  
**Cc:** armand\_scaten@agilent.com  
**Subject:** RE: Software Upgrade  
**Attachments:** Univ of Mass GCMS sw upgrade.pdf

Hi James

As per our earlier correspondence, please see attached your quote for a GCMS software upgrade as well as a Performance Upgrade kit. The performance upgrade will bring your instrument to a LAN interface which is required in order for the instrument to be compatible with this newest revision of software. Also included in with the software upgrade is a PC. If you are submitting this for budgeting I would advise you to add on about \$2,000 if you would like Agilent personnel to complete the installation for you. I also noticed that you have a couple of MS libraries, if you would like to upgrade these as well please let me know and I can add this to your quotation.

As always, please do not hesitate to contact me with any further questions or concerns you may have.

Regards,

**Michelle**

Michelle Tomba | Instrument Online Sales | 302-633-8807

Agilent Technologies | Life Science and Chemical Analysis | [www.agilent.com/chem](http://www.agilent.com/chem)

---

**From:** Hanchett, James (DPH) [<mailto:James.Hanchett@state.ma.us>]  
**Sent:** Monday, July 26, 2010 10:14 AM  
**To:** TOMBA,MICHELLE (A-LittleFalls,ex1)  
**Subject:** Software Upgrade

Hi Michelle,

I hope this information is correct. If not please Email.

1. The model and serial numbers of both the GC and MS.

	Model #	Serial #
GC 8890	G1530A	US00025881
MS 5973	G1098A	US82311418

2. What is your current revision of software as well as your current operating system on your PC?

PC uses Windows NT 4.00.1381

Software	Revision
G1033A	D.01.00
G10391B	B.00.00
G1701BA	B.01.00

3. How is the instrument currently interfaced to the PC?

We use the GPIB interface.

4. Are there any auto samplers on the GC?

	Model #	Serial #
Tower	G2613A	US80204422
Tray	G2614A	US84302348

5. Are you using any headspace samplers, purge and trap, thermal desorption units, etc.?

No.

6. We will need to know the firmware revision of each GC.

Version A.03.05

7. Would you like me to include time for Agilent personnel to perform the install for you?

Not at this time. We are submitting this proposal for funding in the future.

Thanks,

Jim Hanchett  
Amherst Drug Lab  
Room N251 Morrill I  
637 North Pleasant Street  
Amherst, MA 01003